

Lenti-Virus Protocol

张端午

A. For packing virus in 6-well plate: (if in 12-well plate, reduce all to 1/2)

1. Mix the following plasmids in a 1.5 ml eppendorf tube,
1.5 μg Lenti-vector + 1.5 μg packing plasmids = 3 μg total
0.5:0.3:0.2 (pMDL = 0.75 μg , VSV-G = 0.45 μg , REV = 0.3 μg)
add 200 μl OPTI-MEM, mix well.
2. Mix 7.5 μl LF2K in 200 μl OPTI-MEM, incubate for 5 min.
3. Mix the LF2K mixture with plasmids mixture, incubate for 15-20 min.
4. During the incubation time, trypsinize 293T:

Add 1 ml 0.05% Trypsin & EDTA to a full 100 mm plate ($\sim 1.2\text{-}1.4 \times 10^7$ cells),
Add 5 ml fresh medium to terminate reaction,
Count cells and adjust to a concentration of 2×10^6 cells / ml.
5. Aliquot 1 ml cells to a well of 6-well plate (2×10^6 cells / well), then add another 1 ml medium.
6. Add transfection mixture to the well, mix well, and move back to 37C incubator quickly.
7. 8-12 h later, change medium (2 or 2.5 ml), and incubate for another 24-36 h.
8. Harvest virus. (virus can be stored at 4C for one week, otherwise, aliquot virus and store at -80C for further use)

B. Infection:

9. Seed cells that need to be infected in a 12-well or 6-well plate. The best confluence is 40-70%, depending on the different applications.
10. For 12-well plate: (if in 6-well plate, duplicate all)
Infect mouse cells: 400-800 μl virus + 800-400 μl fresh medium = 1.2 ml total
Infect human cells: 100-500 μl virus + 1100 μl -700 μl
11. Add 8-10 $\mu\text{g/ml}$ polybrene, spin for 30 min at 2500 rpm at 37C.
12. 12-24 h later, change medium.